



## IMPORTANT NOTICE FOR USERS OF LEGACY DRUGPAK PRODUCTS

*If you are currently running DrugPak version 5.5 or later, you can disregard this document. Instead see the enclosed "Upgrade Guide" document.*

DrugPak program updates are intended for members of the DrugPak Extended Support plan who remain up-to-date with each major release. It is not possible to skip major releases. If your current DrugPak installation is earlier than version 5.5, then you must install each interim major release and perform the corresponding database upgrade before installing the next release. This document will guide you through that process.

Before installing any program updates, please identify the version of DrugPak you are currently running. Log into DrugPak, and from the main menu, select "Help | About" to see the Version.

### **MAKE A BACKUP OF YOUR CURRENT DRUGPAK PROGRAMS AND DATA BEFORE PROCEEDING!**

*All of the following instructions assume that your CD-ROM drive is drive D:\*

#### **If your current version is DRUG-PAK for DOS, you must upgrade to DrugPak v2.3**

- 1) On the DrugPak 6.0 CD, run "D:\OLD\DP16 Installs\DPSETUP.EXE." Follow the on-screen instructions to install DrugPak v.2.3.
- 2) Run DrugPak v2.3 at least once (login with the username MASTER and no password).
- 3) On the DrugPak 6.0 CD, run "D:\OLD\DP16 Installs\MIGRATE1\MIGRATE.EXE." Follow the on-screen instructions to migrate your data from DRUG-PAK (DOS) to DrugPak v2.3.
- 4) After completing the data migration, proceed to the next section, below.

#### **If your current version is DrugPak v2.x, you must upgrade to DrugPak v3.3**

- 1) On the DrugPak 6.0 CD, browse to the folder "D:\OLD\DP3.2 Installs"
- 2) Run the installation program for each of the modules you currently have installed. DrugPak will install into a new directory (... \Scanlon\dp32)
- 3) Run DrugPak v3.3 at least once and log in. You will be prompted to undertake the database migration.
- 4) After completing the data migration, proceed to the next section, below.

#### **If your current version is DrugPak v3.x, you must upgrade to DrugPak v4.0**

- 1) On the DrugPak 6.0 CD, browse to the folder "D:\OLD\DP4.0 Installs"
- 2) Run the installation program for each of the modules you currently have installed. DrugPak will install itself into the ... \Scanlon\dp32 directory.
- 3) Log into DrugPak v4.0, and follow the on-screen instructions to complete the database upgrade.
- 4) Then run and log into each additional module that you have installed, and follow the on-screen instructions to complete additional database upgrades, if necessary.
- 5) Proceed to the next section, below.

#### **If your current version is DrugPak v4.x, you must upgrade to DrugPak v5.2**

- 1) On the DrugPak 6.0 CD, browse to the folder "D:\OLD\DP5.2 Installs"
- 2) Run the installation program for each of the modules you currently have installed. DrugPak will install itself into a new folder (... \Scanlon\dp50)
- 3) Log into DrugPak v5.2, and follow the on-screen instructions to complete the database migration.

#### **If your current version is DrugPak v5.x, you must upgrade to DrugPak v5.5**

- 1) On the DrugPak 6.0 CD, browse to the folder "D:\OLD\DP5.5 Installs"
- 2) Run the installation program for each of the modules you currently have installed. DrugPak will install itself into a new folder (... \Scanlon\dp50)
- 3) Log into DrugPak v5.5, and follow the on-screen instructions to complete the database migration.

#### **If your current version is DrugPak v5.5, you must upgrade to DrugPak v6.0**

See the enclosed document "Upgrade Guide."